

# Euronda®

It's all about protection



**mammoth**  
DENTAL SUPPLIES



**0330 332 0645**



info@mammothdental.co.uk

## E10 Autoclave



**\*\* 5 YEAR or 4,500 cycle WARRANTY \*\***

- Large 24 litre chamber
- Built in data logger
- Built in Wi-Fi to transfer data to PC
- Holds 5 trays or 6kg load
- Selectable N S & B cycles
- Patented water & hand-piece oil separator
- Accessible clean & dirty water tanks for easy access to clean reducing bio-film build up
- Unique features enhance reliability allowing user to run cycles when the vacuum pump or steam generator have failed

# ONLY 4,995

RRP £7,000.00

### CYCLE TIME'S

**N ( solid non-vacuum cycle)**

**16 minutes**

**B (Solid, lumen & porous vacuum cycle)**

**24 minutes including drying**


**S (Solid & lumen vacuum cycle)**

**21 minutes including drying**



**GREEN Light**

End cycle, load sterile



**YELLOW Light**

End cycle, load wet for prompt use



**RED Light**

Load not sterile



mammothdental.co.uk

## TEAM UP

with colleagues to  
**BUY 5 & SAVE £2,000**

# ONLY 4,595

each

ALL PRICES PLUS CARRIAGE, INSTALLATION & VAT - E&OE 12/18