



**mammoth**  
DENTAL SUPPLIES



**MORITA**



**0330 332 0645**



**info@mammothdental.co.uk**

# The ULTIMATE Dental Unit

## WHY?

### COMFORT

Hydraulic chair with Super sensitive control with variable backrest height and Slow -speed mode ensure optimum comfort for patients

### LIGHT

LED operating light offering shadow free illumination and motion detector.

### FOOTPRINT

Small footprint size allows for maximum legroom & precise focusing

### INSTRUMENT SETUP

Innovative instrument design allows for single or four handed dentistry, together with anti suck back mechanism to prevent contamination

### CUSDOR

Swivels to ensure easy use together with automatic optical cup refill

### TEMPERATURE CONTROL

A warm water boiler for instrument spray and mouthwash water

### DISINFECTION

Integrated hose & disinfection system ensures highest levels of hygiene

### FILTER SYSTEM

A Micro filter system for fresh water and compressed-air supplies ensures clean air and water without additional filters

### TOUCH CONTROL

LCD non-reflective touch-screen display allows for easy use and function

### ELECTRIC MOTOR

Torx micro motor with LED fibre optics with speed control on display

### FUTURE

USB-connection for camera or Multimedia devices

# #SOARIC

**From**  
**£29,995**

**100+** Years  
A Century of Innovation



RRP £44,250.00



**mammothdental.co.uk**



**5 YEAR WARRANTY**

ALL PRICES PLUS CARRIAGE & VAT - E&OE 12/18

**EXTRA'S:-**

~ KNEE-BREAK CHAIR ~ BUILT IN IMPLANT MOTOR ~  
~ APEX LOCATOR ~ MICROSCOPE ~